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Resdevmed Lda.
Estrada da Marinha, 542
3880-123 Ovar, Portugal

tel (351)256782047
www.boneeasy.com

Custom-made craniomaxillofacial implant in PEEK

The use of the custom-made **iCMF PEEK** medical device should only be carried out by a qualified clinician, with in-depth knowledge and complete command of the specific surgical techniques of implant procedures and craniomaxillofacial prostheses. Before using a BoneEasy (Resdevmed Lda.) brand product, the responsible surgeon/doctor must carefully analyze all information provided by the manufacturer, including indications,

contraindications, warnings, instructions and other information, and respect them. Resdevmed Lda. recommends participation in practical training courses on the use of BoneEasy brand products. Detailed instructions in addition to those contained in this document can be obtained by contacting the manufacturer or its representative.

iCMF PEEK should be used at the discretion of the clinician, who has the obligation to determine whether the product is suitable for the patient and evaluate all relevant circumstances. The clinician is responsible for any direct or indirect complications or harmful situations that may result from an incorrect indication or surgical technique, inappropriate use of the material, overload, lack of asepsis, or failure to observe the explicit safety instructions in the instructions for use. The manufacturer or representative of RESDEVMED Lda. cannot be held responsible for complications associated with use by the doctor, as described above, or by the patient, including

the patient's anatomy and general habits. Resdevmed Lda. declines any responsibility, express or implied. It is also the responsibility of the surgeon/doctor to adequately inform the patient about the function and necessary care, as well as the known risks associated with the product in question.

Description

iCMF PEEK is a custom-made craniomaxillofacial implant made of Poly-ether-ether-ketone. It is developed based on files from a computerized tomography, using computer-aided segmentation and design software and manufactured by Fused Deposition modeling (FDM). The **iCMF PEEK** medical device must be stabilized with fixing screws, not included with this device.

Indications

iCMF PEEK is a custom-made, long-term medical device for use in rehabilitation and/or replacement procedures in craniomaxillofacial areas.

Contraindications

iCMF PEEK should only be used if the volume of cortical bone allows stabilization of the implant. Poor bone quality, difficulty anchoring screws, systemic diseases (diabetes, etc.), reduced immunity, alcoholism, drug addiction and

psychological instability are factors that can contribute to the lack of integration and/or subsequent failure of the implant. Incorrect use of these implants or poor fixation can lead to screw relaxation, fractures and/or implant failure. Exposure to radiation and chemotherapy can affect the functionality of the implant.

Warnings

Incorrect surgical techniques can result in bone loss, patient harm, pain, and partial or complete failure of the medical device.

Steroid or anticoagulant treatments can affect the surgical site and impact the patient's ability to integrate.

Prolonged exposure to or use of bisphosphonate-based medications, especially with chemotherapy, can have a negative impact on implant functionality. A detailed study of the patient's history is recommended, including consultation with the attending physician, before opting for any solution proposed by Resdevmed Lda.

Mobility, bone loss or infection may be indications of implant failure. The decision on the best course of action must be carefully evaluated by the doctor. If implant removal is necessary, closure by first intention should be sought.

Precautions

Proper clinical case planning is crucial to the long-term success of the medical device. During surgery, asepsis rules must be respected. Direct handling of the medical device should be avoided.

iCMF PEEK cannot be used in load situations.

There cannot be micromovements of these devices as they lead to bone resorption, therefore adequate fixation with osteosynthesis screws is essential.

iCMF PEEK is a medical device custom-made for a particular patient, and should therefore not be used on any patient other than the one on which its manufacture is based.

iCMF PEEK was developed for single use. It must not be reused, reprocessed or resterilized. Failure to follow these instructions may compromise the structural integrity of the device and/or lead to its failure, with consequent damage to the patient.

The patient's clinical situation must be carefully

monitored. Do not use the medical device if its original packaging is open, damaged or shows signs of deterioration. See also Contraindications.

Requirements

iCMF PEEK is a medical device designed on a tomography where three-dimensional volume reconstruction is performed. For a good adaptation of

the device, the quality of the tomography is essential. To do this, the presence of metals that can cause effects and that do not allow the volume to be reconstructed correctly must be checked. The

removal of metallic elements that interfere with the three-dimensional reconstruction may be necessary prior to carrying out the examination. The request to the radiology center must include frequent cuts (around 2 mm) unless the volume is very large and causes harmful exposure to the patient, this being a clinical criterion that must be evaluated by the doctor.

Recommendations

During the first week, at least one visit for patient monitoring and professional prophylaxis is recommended. X-rays can be taken four months after surgery to assess the condition of the tissues and the medical device, unless complications arising from implantation require early assessment tests.

Possible adverse effects

Complications that may occur from the use of the implant include (but are not limited to): pain, discomfort, delayed healing, paresthesia, hyperesthesia, edema, hemorrhage, hematoma, infection, inflammation, need for additional surgery or removal, lack of integration, mechanical failure, exposure, bone loss and implant loss. Other adverse effects may also occur as a result of iatrogenic factors or patient response.

Implant removal should be considered whenever the location where it was placed shows signs of being compromised in such a way that it cannot be controlled by postoperative treatments.

Report any adverse situation recorded and not presented in this document to the manufacturer.

Technical information

iCMF PEEK is an implantable medical device that requires proper planning. Incisions and detachments that allow intimate contact with the bone must be previously planned, it must be assessed whether the added volume is sufficient to be contained by the soft tissues and complementary surgical techniques that may be necessary to complement the procedure must be foreseen, such as such as microvascularized flaps, grafts or the previous use of tissue distension balloons. Excessive torques on the screws must be avoided, either because they can create opposite vector forces on the screws, leading to their loosening, or due to the mechanical limits of the implant, which can lead to its fracture.

Sterilization

iCMF PEEK is sterilized by ethylene oxide. The packaging will serve as a sterilization barrier until the expiration date indicated on the box. This medical device was developed for single use and should not be re-sterilized.

Single use

iCMF PEEK must not be reused. Reusing a single-use medical device that has been in contact with blood, bone, tissue, body fluids or other contaminants may lead to harm to the user. Possible risks associated with reusing a single-use device include, but are not limited to, mechanical failure and transmission of infectious agents. **iCMF PEEK** is a medical device custom-made for a particular patient, and should therefore not be used on any patient other than the one on which its manufacture is based.

Packaging

iCMF PEEK was cleaned and packaged in a controlled environment. It is supplied in a multipack. The outer label contains information about the batch number that must be recorded in the patient's clinical record to guarantee complete traceability of the product. The manufacturer also provides extra labels on the packaging that can be placed on the clinical record for the same purpose.

Storage

iCMF PEEK must be stored at room temperature and protected from external damage. If contaminated with blood or other fluids, it must be disposed of in accordance with rules for biological contamination waste.

Patient information

It is the responsibility of the surgeon/doctor to adequately inform the patient about the function and necessary care, as well as the known risks associated with the product in question.

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