IFU - Mesh4U

Instruction for use - Mesh4U

RG.PR.20.04

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www.boneeasy.com

Mesh4U

Custom-made titanium mesh for guided bone regeneration

Mesh4U is a custom-made medical device, manufactured specifically under medical prescription and under the responsibility of the prescriber, with specific design features intended for use exclusively in a particular patient. **Mesh4U** should only be used by a qualified clinician with thorough knowledge and mastery of the specific surgical techniques of guided bone regeneration procedures. Before using a BoneEasy (Resdevmed Lda.) branded product, the responsible surgeon/physician should carefully review all information provided by the manufacturer, including indications, contraindications, warnings, instructions and other relevant information. Detailed instructions other than those contained herein may be obtained by contacting the manufacturer or its representative.

Mesh4U should be used at the discretion of the clinician, who has an obligation to determine if the product is suitable for the patient and to evaluate all relevant circumstances. The clinician is responsible for any direct and/or indirect complications or detrimental situations that may result from erroneous indication or surgical technique, misuse of material, overload, lack of asepsis, or failure to follow explicit safety instructions in the instructions for use. The manufacturer or representative of Resdevmed Lda. cannot be held responsible for complications associated with use by the physician as described above or to the patient, including the patient's anatomy and general habits. Resdevmed Lda. disclaims any liability, express or implied. It is also the surgeon/doctor's responsibility to adequately inform the patient about the function and necessary care, as well as the known risks associated with the product.

Description

Mesh4U is a custom-made metallic barrier mesh made of medical grade titanium (Ti6Al4V). It is developed based on files from computed tomography, using computer aided design and segmentation software, and manufactured by laser sintering (SLM Selective Laser Melting).

Mesh4U must be stabilized with screws not included with this device.

Indications

Mesh4U is a long-term implantable medical device suitable for guided bone regeneration surgical procedures.

Contraindications

The clinician should be careful in selecting the patients indicated for guided bone regeneration procedures, in selecting the most appropriate postoperative treatments and should also be critical during the **Mesh4U** design process. **Mesh4U** is not designed to handle high loads. For this reason the load on the individualized mesh should be zero, with removable prostheses resting on the tissues that cover it to avoid its collapse and not to compromise the healing process. **Mesh4U** can be designed with pillars for provisionalization, for aesthetic use only and not for loading. This medical device was not designed for joints (eg, mandibular temporal joint reconstruction). The long-term effectiveness and safety of combining **Mesh4U** with bone fillers is not yet established, so when performed in combination, the clinician should follow individual material indications.

The following conditions may affect the medical device performance: deficient oral hygiene, organ transplants, cardiovascular disease, hypertension, thyroid or parathyroid disease, alcohol or drugs, tobacco, anticoagulant, anticonvulsant, immunosuppressive therapies, with inflammatory and/or infectious active processes of the oral cavity, cardiovascular failure, coronary disorders, arrhythmia, chronic pulmonary or respiratory diseases, gastrointestinal disorders, hepatitis, intestinal inflammation, chronic renal failure and urinary system disorders, endocrine disorders, haematological problems, anemia, leukemia, clotting problems, osteoporosis or musculoskeletal arthritis, stroke, neurological disorders, mental retardation or paralysis, chemotherapy, bisphosphonate administration (especially orally and intravenously), presence of periodontal disease, titanium allergy, chronic acute infectious diseases, chronic subacute maxillary osteitis, systemic diseases, endocrine disorders, pregnancy, lactation, renal failure, fibrous dysplasia, hemophilia, neutropenia, steroid use and diabetic disease.

Other normal contraindications common to all oral surgery interventions should also be considered.

Warnings

Inadequate surgical techniques can result in bone loss, harm to the patient, pain, and partial or complete medical device failure.

A detailed study of the patient's history, including consultation with the attending physician, is recommended before opting for any solution available at Resdevmed Lda.

The clinical condition of all patients should be continuously monitored and, if necessary, the medical device removed.

When combining **Mesh4U** with any type of bone defect replacement, all indications for each material should be considered individually.

Mesh4U can be used in conjunction with endosseous implants, under prescriber doctor indications.

Due to the conductivity of the metal, electrosurgery can result in tissue damage, patient injury and implant failure and should be used with caution.

See also Contraindications.

Precautions

Proper case planning is crucial for the long-term success of the medical device.

During surgery aseptic rules must be respected. Direct handling of the medical device should be avoided.

Mesh4U should not be used under load situations.

Mesh4U is a custom-made medical device for a particular patient, so it should not be used on a patient other than the one for whom it is manufactured.

Mesh4U is designed for single use only. It must not be reused, reprocessed or resterilised. Failure to follow these instructions can compromise the structural integrity of the device and/or lead to device failure, with consequent patient damage.

The clinical situation of the patient should be carefully monitored.

Do not use the medical device if its original packaging is open, damaged or showing signs of deterioration.

See also Contraindications.

Recommendations

An oral hygiene plan should be prescribed by the clinician, which may include mechanical and chemical control of plaque and instructions for oral hygiene.

Antibiotic therapy is recommended at the discretion of the clinician.

During the first week, at least one visit is recommended for patient monitoring and prophylaxis.

X-rays may be taken after surgery to assess tissue and medical device status, unless complications from implantation require early screening.

Medical device removal should be considered in case of exposure, complications that cannot be controlled by standard postoperative treatments, tissue inflammation or evidence of infection, but always at the discretion of the clinician.

Possible adverse effects

Complications that may occur from the use of this medical device include (but are not limited to): pain, discomfort, edema, bruising, inflammation, thermal sensitivity, infection, exfoliation, perforation or abscess formation, hyperplasia, gingival irregularities, complications associated with anaesthesia, mechanical failure of the medical device or exposure. Other adverse effects may also occur as a result of iatrogenic factors or patient response.

Removal of the mesh should be considered whenever the site where it is placed shows signs of compromise in such a way that it cannot be controlled by postoperative treatments.

Report to the manufacturer any adverse events recorded and not shown in this document.

Technical information

Mesh4U is an implantable medical device that requires proper planning.

The following considerations are proposed by Resdevmed Lda. However, it is important to remember that **Mesh4U** implantation should only be performed by qualified clinicians, with thorough knowledge of the specific surgical techniques of guided bone regeneration procedures.

For **Mesh4U** placement during surgery:

1. Keep the sterile field throughout the procedure.
2. Minimize saliva or any other source of contamination to material and surgical site.
3. Gently open the outer blister and remove the inner blister containing the sterile **Mesh4U** into the sterile field. Carefully remove the medical device from the internal blister.
4. Place the medical device in the area to be treated.
5. The stability of the medical device is ensured by fixing with suitable screws.
6. To aid in the regenerative process, the medical device should be combined with bone substitutes and covered with a resorbable membrane.

Removal of the medical device is the clinician's discretion: depending on the type of application, different time windows may be recommended for removal of implanted materials: three to nine months or until bone regeneration is complete to place the implants; four to twelve weeks for transgingival healing.

Sterilization

Mesh4U is sterilized by ethylene oxide. The packaging will serve as a sterilization barrier until the expiration date indicated on the box. This medical device is designed for single use only and should not be resterilized.

Single use

Mesh4U should not be reused. Reuse of a single use medical device that has been in contact with blood, bone, tissues, body fluids or other contaminants may lead to damage to the user. Possible risks associated with reusing a single use device include, but are not limited to, mechanical failure and transmission of infectious agents. **Mesh4U** is a custom-made medical device for a particular patient, so it should not be used on a patient other than the one on which it is manufactured.

Packaging

Mesh4U has been cleaned and packed in a controlled environment. It is supplied in multiple packaging. The outer label contains information about the batch number that should be recorded on the patient's clinical record to ensure complete traceability of the product. The manufacturer provides extra labels, available on the package, which can be placed on the medical record for the same purpose. One of the extra labels provided must be given to the patient.

Storage

Mesh4U should be stored at room temperature and protected from external damage.

Disposal, in the case of post-surgical removal, must follow rules for disposal of contaminants with blood.

Disposal of produced parts, without contact with biological contamination, must follow titanium disposal rules.

Patient Information

It is the surgeon/physician's responsibility to adequately inform the patient of the necessary function and care as well as the known risks associated with the product.

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